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## **Instructions for Use of TPP Vacuum Filtration «rapid»-Filtermax**

*The TPP filtration system «rapid»-Filtermax is suitable for sterile filtration of cell culture media, sera and aqueous solution.*

*The 0.22 µm polyether sulfone (PES) membrane ensures high flow rates with low protein binding and foam formation. The «rapid»-Filtermax is available as a filtration unit (bottle top filter incl. reservoir flask) as well as a separate bottle top filter or flask, with a DIN GL45 thread.*

*The filtration system is available in volume sizes 150, 250, 500 and 1000 mL. Each system component is identical in construction, available separately, and can be combined with each other.*

*The square foot of the reservoir ensures a safe stand during sterile filtration.*

*The «rapid»-Filtermax is for single use only.*

### **Note:**

*Follow the national regulations for handling biological materials and wear the appropriate protective clothing.*

*During the work process, be sure to observe the rules of aseptic technique.*

### **Instructions:**

- *Inspect packaging and product for integrity. Use only materials in perfect condition.*
- *Check the expiry date (EXP) on the label and packaging. Use only products with a valid EXP.*

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### **Advice and Recommendations:**

- *The volume scale is an indication of the filling volume, not an absolute value.*
- *Alcohol 90% with mechanical stress (rubbing) may dissolve or blur the white color of the print.*

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## Information

- *The user should provide the following materials: Vacuum source; vacuum tubing.*
- *Centrifuge or prefilter solutions with high particle densities as a priority to prevent filter clogging.*
- *Attach solutions with a low particle density and a high filtration throughput to another reservoir for additional filtration. Please observe the rules of aseptic technique.*
- *Once the pouch is open, keep the product in a sterile area to ensure sterility.*
- *The flask is not suitable for freezing.*
- *Do not autoclave product.*
- *This product is not for direct patient care or diagnostic procedures, but for research use only.*

## To avoid possible injury due to implosions:

- *Do not exceed 700 mm Hg differential vacuum at 25 °C*
- *Bottle top filter: Use only sterile plastic or glass bottles that are approved for sterile filtration under vacuum. These should be free of any damage (risk of implosion). Observe the maximum total pressure.*
- *Wear eye-protection and perform the filtration in special safety units.*
- *Do not use hand-held vessels while filtering under vacuum.*

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**Technical Data:**

**Materials**

<i>Filter-membrane</i>	<i>Polyethersulfone (PES), type TPP fast flow</i>
<i>Flask and top</i>	<i>PS</i>
<i>Screw cap</i>	<i>PE</i>
<i>Vacuum connector</i>	<i>PP</i>

<b>Technische Daten</b>	<b>Technical Data</b>	<b>Werte / Values</b>
Max. Gebrauchs-Temperatur °C	<i>Max temp in use °C</i>	max. 45
Vakuum mm Hg bei 25 °C	<i>Vacuum mm Hg at 25 °C</i>	700
Totvolumen ml	<i>Dead volume mL</i>	2.3
Vakuumschluss mm	<i>Vacuum connection mm</i>	6 – 10
Wasserflussrate bei 25 °C	<i>Water flow rate at 25 °C</i>	22 ml/min/cm <sup>2</sup> 0.07 mPa (0.7 bar)
Proteinbindung µg/cm <sup>2</sup>	<i>Protein binding µg/cm<sup>2</sup></i>	< 20

		<b>System</b> <i>System</i>	<b>Flaschenaufsatzfilter</b> <i>Bottle Top Filter</i>	<b>Flasche</b> <i>Flask</i>
<b>Produkt Nr.</b>	<b>Product No.</b>	<b>99150</b>	<b>99155</b>	<b>99157</b>
<b>Abmessungen</b>	<b>Measurements</b>			
Höhe mm	<i>Height mm</i>	103	57	63
Breite x Länge mm	<i>Width x Length mm</i>	93 x 93	90 x 89	93 x 93
Filtergrösse cm <sup>2</sup>	<i>Filter size cm<sup>2</sup></i>	49	49	---
Porengrösse µm	<i>Pore size µm</i>	0.22	0.22	---
Volumen ml	<i>Volume mL</i>	150	150	150

		<b>99250</b>	<b>99255</b>	<b>99257</b>
<b>Produkt Nr.</b>	<b>Product No.</b>			
<b>Abmessungen</b>	<b>Measurements</b>			
Höhe mm	<i>Height mm</i>	143	75	87
Breite x Länge mm	<i>Width x Length mm</i>	93 x 93	90 x 89	93 x 93
Filtergrösse cm <sup>2</sup>	<i>Filter size cm<sup>2</sup></i>	49	49	---
Porengrösse µm	<i>Pore size µm</i>	0.22	0.22	---
Volumen ml	<i>Volume mL</i>	250	250	250

		<b>99500</b>	<b>99505</b>	<b>99507</b>
<b>Produkt Nr.</b>	<b>Product No.</b>			
<b>Abmessungen</b>	<b>Measurements</b>			
Höhe mm	<i>Height mm</i>	213	111	121
Breite x Länge mm	<i>Width x Length mm</i>	93 x 93	90 x 89	93 x 93
Filtergrösse cm <sup>2</sup>	<i>Filter size cm<sup>2</sup></i>	49	49	---
Porengrösse µm	<i>Pore size µm</i>	0.22	0.22	---
Volumen ml	<i>Volume mL</i>	500	500	500

		<b>99950</b>	<b>99955</b>	<b>99957</b>
<b>Produkt Nr.</b>	<b>Product No.</b>			
<b>Abmessungen</b>	<b>Measurements</b>			
Höhe mm	<i>Height mm</i>	285	143	160
Breite x Länge mm	<i>Width x Length mm</i>	111 x 111	108 x 108	111 x 111
Filtergrösse cm <sup>2</sup>	<i>Filter size cm<sup>2</sup></i>	69	69	---
Porengrösse µm	<i>Pore size µm</i>	0.22	0.22	---
Volumen ml	<i>Volume mL</i>	1000	1000	1000

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**Additional:**

*Instructions for use, chemical resistance lists, and quality certificates for each product are available for download from [www.tpp.ch](http://www.tpp.ch).*