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Instruction for Use TPP TubeSpin® Bioreactor with Septum

The TubeSpin® Bioreactor with Septum combines the features of the standard TubeSpin® Bioreactor #87050 with the requirements of laboratory automation. It enables multiple access points for injection needles and pipette tips within laboratory automation workflows, facilitating the exchange of sterile media, solutions, or the sampling of cell-based products.

The TubeSpin® bioreactor with a silicone septum is available with a filter screw cap with openings (Figure 1, #1) above a PTFE filter membrane, and with or without barcode 128.

The TubeSpin® bioreactor with septum is for single use only.

Safety instructions:

1. Handling and Safety

When dealing with biological material, adherence to national regulations is mandatory. Furthermore, it is essential to ensure compliance with the established biological safety level of the laboratory, along with the safety data sheets and the manufacturer's instructions for use.

Wear your personal protective equipment.

2. Risk of contamination

Always follow the rules for aseptic techniques and good laboratory practices. Open the packaging only immediately before use and use only items that are visually flawless. Dispose of any products that show visible defects or signs of contamination.

3. Storage

Store the TPP products at room temperature and protect them from direct UV radiation.

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Instructions:

- *Check the expiration date (EXP) on the label and packaging. Only use products with a valid EXP date.*
- *Before use, verify that the packaging is intact, as the consumable is only considered sterile if the packaging is undamaged.*
- *Open the packaging in a sterile environment and remove the product.*
- *Open the bioreactor or use the access through the septum and fill the tube with medium and inoculum according to standard laboratory procedures. Ensure the optimal filling volume.*
- *Close the cap and ensure it is properly seated on the tube to minimize the risk of contamination.*
- *Continuous gas exchange takes place through the openings above the hydrophobic 0.22 µm membrane. If the PTFE filter membrane gets wet, the gas exchange is temporarily reduced.*
- *The silicone septum with cross slit is accessible for needles and tips with a diameter of <2.1 mm, fill according to normal laboratory routine.*
- *Multiple access through the septum with needles and tips are possible keeping the compliance of the aseptic work. No maximum values or numbers were determined for multiple access.*
- **Laboratory automation:** *Refer to the operating instructions of the manufacturer. Observe the optimum filling quantity (see technical data).*
- *Barcode 128 on each tube for efficient automated tracking.*
- **Cultivation:** *Shake the cells in appropriate incubation shakers. Observe the cultivation conditions of the cells (temperature, CO₂ and humidity control, shaking frequency).*

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
Figure 1

Remarks and Recommendations:

- The graduation serves as a reference for the fill volume and is not an absolute value.
- Alcohol (90%) combined with mechanical stress (rubbing) may dissolve or blur the white labeling.

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Usage and storage at sub-zero temperatures

- **TPP Note:** Storage at subzero temperatures is not an intended application for this TPP product. TPP recommends transferring samples for these conditions into TPP cryo tubes, which are specifically designed for this purpose.
- Do not expose the consumables to liquid nitrogen. They may become damaged or burst when thawed.
- When using bioreactors made of polypropylene (PP) at subzero temperatures, a preliminary application test should be conducted to ensure that the material withstands the specific requirements.

Additional:

Instructions for use, chemical resistance lists, and quality certificates for each product are available for download from www.tpp.ch.

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Disclaimer:

TPP products are intended solely for use in research and are not suitable for clinical, diagnostic, or therapeutic applications.

TPP assumes no liability for damages or defects resulting from improper or unauthorized use. Users are responsible for storing, handling, and using the products in accordance with the provided instructions.

TPP's recommendations are general guidelines and may not cover all possible scenarios. TPP is not liable for indirect, incidental, consequential, or special damages arising from the use or misuse of these products.

Swiss law governs these terms of use and any resulting legal matters.