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## Instruction for Use TPP Cryo Tube

*TPP cryo tubes are made of polypropylene (PP) with sealing lip in the cap, external thread tube and star feet are designed for storage of biological material and samples down to -196 °C. The screw caps with external thread reduce the risk of contamination and cross-contamination.*

*The cryo tube is for single use only.*

### **Safety instructions:**

#### 1. Handling and Safety

*When dealing with biological material, adherence to national regulations is mandatory. Furthermore, it is essential to ensure compliance with the established biological safety level of the laboratory, along with the safety data sheets and the manufacturer's instructions for use.*

*Wear your personal protective equipment.*

#### 2. Risk of contamination

*Always follow the rules for aseptic techniques and good laboratory practices. Open the packaging only immediately before use and use only items that are visually flawless. Dispose of any products that show visible defects or signs of contamination.*

#### 3. Storage

*Store the TPP products at room temperature and protect them from direct UV radiation.*

### **Notes on storing in nitrogen**

*TPP recommends storing the tubes exclusively in the gas phase of liquid nitrogen (LN<sub>2</sub>). Storage in the liquid phase could allow LN<sub>2</sub> to enter the tubes through small openings, such as for example improperly closed caps. During the thawing process, rapid evaporation will generate internal high pressure, potentially causing the tubes to explode and release biologically hazardous substances.*

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### Instruction

- *Check the expiration date (EXP) on the label and packaging. Only use products with a valid EXP date.*
- *Before use, verify that the packaging is intact, as the consumable is only considered sterile if the packaging is undamaged.*
- *Open the packaging in a sterile environment and remove the product.*
- *For a firm and secure stand while working, place the cryo tubes in the cryo tube rack (# 99016).*
- *Unscrew the lid of the cryogenic tube and fill the tube according to the usual laboratory routine.*
- *Please note the maximum filling quantity of 90% <sup>[1]</sup> to allow for expansion of the sample volume (negative thermal expansion of water)*
- *Close the cap tightly. The thread must be dry to close the tube, moisture will affect the tightness. Avoid overtightening the cap and damaging the sealing lip.*
- *Use only cryo boxes (or freezers) designed for freezing procedures. Unsuitable boxes and racks have an insulating effect and may cause uneven freezing, resulting in loss of cell viability and/or damage to the tube.*
- **Centrifugation:** *Use suitable rotors or centrifuge adapters and always follow the safety instructions of the centrifuge manufacturer*
- **Factors affecting the mechanical load capacity of tubes:**
  - *Shape and material of the tubes*
  - *Fit in the centrifuge adapter*
  - *Temperature, centrifugation time, and relative centrifugal force (RCF)*
  - *Chemical and physical properties of the centrifuged material*
  - *Type of rotor (fixed-angle or swing-out rotor)*

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- *Before regular use of the product, perform a test run with the chosen settings to verify its suitability.*

#### **Storage in Minus Temperatures:**

- *When storing the cryo tubes in the subzero range, first test the many influencing factors under routine conditions.*
- *Temperatures below 0 °C result in reduced elasticity for polypropylene. Do not subject frozen tubes to mechanical stress, such as dropping or hitting it.*
- *Use only "freezing containers" for controlled freezing of cells before transferring the cryo tubes to the gas phase of a nitrogen tank.*

#### **Advice and Recommendations:**

- *The volume scale is an indication of the filling volume, not an absolute value.*
- *A labeling area and colored cap inserts allow for sample identification.*
- *An example of a cell freezing protocol can be found on our website at [www.tpp.ch](http://www.tpp.ch).*

#### **Additional:**

*Instructions for use, chemical resistance lists, and quality certificates for each product are available for download from [www.tpp.ch](http://www.tpp.ch).*

#### **Accessories:**

- *Cryo Box L, # 99014 / Cryo Box S, # 99015*
- *Cryo Rack # 99016: The star foot in combination with the TPP rack allows one-handed operation.*
- *Colored cap inserts for identification*

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**Instructions:**

Designed for the storage of cryo tubes and most other similarly designed vials at ultra-low temperatures.

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### **Disclaimer:**

*TPP products are intended solely for use in research and are not suitable for clinical, diagnostic, or therapeutic applications.*

*TPP assumes no liability for damages or defects resulting from improper or unauthorized use. Users are responsible for storing, handling, and using the products in accordance with the provided instructions.*

*TPP's recommendations are general guidelines and may not cover all possible scenarios. TPP is not liable for indirect, incidental, consequential, or special damages arising from the use or misuse of these products.*

*Swiss law governs these terms of use and any resulting legal matters.*

---

[1] Negative thermal expansion of water ~ 10%