## **Gebrauchsanleitung TPP PCV-Röhrchen (Packed Cell Volume)**

### **TPP «easy read» Messinstrument**

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## **Instruction for Use TPP PCV (Packed Cell Volume) tube TPP «easy read» Measuring Device**

*The PCV tube, together with the "easy read" measuring device, allows easy and fast measurement of biomass in cell culture, microbiology and suspension cultures of fungi and yeasts. You need a centrifuge with swing-out rotor.*

*The PCV tube is for single use only.  
The "easy read" measuring device is reusable.*

#### **Safety instructions:**

##### **1. Handling and Safety**

*When dealing with biological material, adherence to national regulations is mandatory. Furthermore, it is essential to ensure compliance with the established biological safety level of the laboratory, along with the safety data sheets and the manufacturer's instructions for use.*

*Wear your personal protective equipment.*

##### **2. Risk of contamination**

*Always follow the rules for aseptic techniques and good laboratory practices. Open the packaging only immediately before use and use only items that are visually flawless. Dispose of any products that show visible defects or signs of contamination.*

##### **3. Storage**

*Store the TPP products at room temperature and protect them from direct UV radiation.*

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
### Instructions TPP PCV tube

- Check the expiry date (EXP) on the label and packaging. Use only products with a valid EXP.
- Open the packaging and remove the product.
- Transfer a sample of 100 - 1000  $\mu$ l of a mixed cell suspension into the PCV cell measuring tube.
- Close the tube with the cap # 87008 if necessary.
- Centrifuge the tube for 1 minute at 2'500 x g in a swing-out rotor for 1.5 | 2 mL vessels.
- Read the volume of the cell pellet as described in the TPP "easy read" measuring instrument instructions.
- Discard the PCV tube.

### Instruction "easy read" Measuring Device

- After centrifugation, insert the PCV tube ① into the opening provided in the "easy read" after centrifugation.
- Move the magnifying glass ③ until the height of the cell pellet in the capillary is equal to the cross.
- Read the cell pellet volume in  $\mu$ l with the "easy read" measuring device.

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### **Calculation PCV (%) value:**

- The measured volume is expressed in % of the volume of the transferred suspension sample. It corresponds to the total biomass of the culture.

$$\text{PCV}(\%) = \frac{\text{volume cellpellett}}{\text{volume probe}} \times 100$$

### **Cleaning and disinfection:**

#### **“easy read” measuring device: Body:**

- Wipe with a soft cloth dampened with a small amount of neutral detergent to remove stains.
- Then wipe the “easy read” clean with a cloth dampened with lukewarm water. Be careful not to touch the magnifier when cleaning.
- Avoid abrasive cleaning agents, scratching metal brushes, agents containing halogens or chlorides: risk of corrosion.
- For disinfection, TPP recommends 70% ethanol.

#### **“easy read” measuring device: magnifying glass**

- Use compressed air to remove dirt and dust adhering to the surface.
- Then take a piece of lens paper and fold it around your finger to form a tip. Be sure to use lens paper for this step, not lab tissues or paper towels. Traditional tissues contain loose, coarse fibers that can scratch or leave marks on the lens.
- Put a small amount of 70% ethanol on the tip of the lens paper.
- When cleaning the lens surface, wipe in a circular motion from the inside to the outside.

Allow the “easy read” to dry completely after cleaning/disinfection. TPP recommends dust-free and dry storage of the device.

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**Comments and recommendations:**

- *Sample dilution is not required even at high cell densities. A reduction of the selected sample volume is sufficient.*
- *The PCV method does not provide direct information on the cell viability of the culture.*
- *Cell pellet cannot be removed from capillary*

**Additional:**

*Instructions for use, chemical resistance lists, and quality certificates for each product are available for download from [www.tpp.ch](http://www.tpp.ch).*

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**Technical Data:**

**Materials**

Measuring device aluminum  
 Magnifying lens glass

Abmessungen	Measurements	87010
Volumengraduierung µl	Volume graduation µL	0.0 – 5.0
Länge x Breite x Höhe mm	Length x width x height mm	320 x 15 x 13.5
Lupe:	Magnifying lens:	
Länge x Breite x Höhe mm	Length x width x height mm	30 x 33 x 38

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### **Disclaimer:**

*TPP products are intended solely for use in research and are not suitable for clinical, diagnostic, or therapeutic applications.*

*TPP assumes no liability for damages or defects resulting from improper or unauthorized use. Users are responsible for storing, handling, and using the products in accordance with the provided instructions.*

*TPP's recommendations are general guidelines and may not cover all possible scenarios. TPP is not liable for indirect, incidental, consequential, or special damages arising from the use or misuse of these products.*

*Swiss law governs these terms of use and any resulting legal matters.*