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Instructions for Use TPP Tissue Culture Flasks

TPP tissue culture flasks are intended for tissue culture in manual use. For optimal cell adhesion and growth, only the bottom of the tissue culture flask is opto-mechanically activated. The angled neck reduces the risk of medium contacting the inside of the screw cap by spillage and allows excellent access for cell scrapers and serological pipettes.

TPP tissue culture flasks are available with filter or VENT screw caps.

The tissue culture flask is for single use only.

Safety instructions:

1. Handling and Safety

When dealing with biological material, adherence to national regulations is mandatory. Furthermore, it is essential to ensure compliance with the established biological safety level of the laboratory, along with the safety data sheets and the manufacturer's instructions for use.

Wear your personal protective equipment.

2. Risk of contamination

Always follow the rules for aseptic techniques and good laboratory practices. Open the packaging only immediately before use and use only items that are visually flawless. Dispose of any products that show visible defects or signs of contamination.

3. Storage

Store the TPP products at room temperature and protect them from direct UV radiation.

Instruction

- *Check the expiration date (EXP) on the label and packaging. Only use products with a valid EXP date.*

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Filter Screw Cap:

- *Close the filter screw cap with ventilation holes. Continuous gas exchange occurs through the integrated 0.22 µm hydrophobic membrane (Figure 1). If the membrane becomes wet, gas exchange is temporarily reduced.*


Figure 1: Filter-screw-cap closed tightly

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
VENT Screw Cap:

- For a gas exchange, turn the VENT cap into the ventilation position at 12 o'clock. You will hear/notice a click.


- Visual check:
Rectangle at 12 o'clock (Position 1)

During incubation of cells, the rectangle must be in position 1 / 12 o'clock.

- To stop gas exchange, rotate the cap 1/4 turn clockwise beyond the vent position. The flask is now tightly sealed and there is no gas exchange.

- Visual check:
Rectangle at 3 o'clock. (Position 2)

- For a transfer of the flask within in the laboratory (laminar flow ↔ incubator), TPP recommends closing the VENT screw cap tightly.

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Optimization of Cell Growth:

- To achieve even cell growth on the surface, please consider the following points:
 - Thoroughly resuspend the cells in the medium to prevent cell clumping.
 - Avoid foam formation during resuspension and seeding.

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Media volume in adherent cell cultures:

- *Use 0.2 - 0.5 mL per cm² of growth area, which corresponds to a medium height of 2 - 5 mm [1].*

The medium height - and thus the volume in the culture vessel - is a critical factor for the oxygen supply to the cells and influences the oxygen transfer rate, OTR (Gstraunthaler et al. 1999).

- *Before regular use of the product, perform a test run with the chosen settings to verify its suitability.*

Advice and Recommendations:

- *For labelling there are bilateral inscription fields*
- *The volume scale is an indication of the filling volume, not an absolute value.*
- *Alcohol (90%) combined with mechanical stress (rubbing) may dissolve or blur the white labeling.*
- *The stacking rim allows a stable stacking of several flasks on top of each other, also in combination with other TPP flasks of the same size.*
- *For homogeneous cultivation of cells, it is essential to level the incubator and provide a continuous supply of medium.*

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- *Avoid and/or eliminate sources of vibration near the incubator to ensure uniform cell growth.*

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Usage and storage at sub-zero temperatures

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- **TPP Note:** Storage at subzero temperatures is not an intended application for this TPP product. TPP recommends transferring samples for these conditions into TPP cryo tubes, which are specifically designed for this purpose.
- Do not expose the consumables to liquid nitrogen. They may become damaged or burst when thawed.
- When using tissue culture flasks at subzero temperatures, a preliminary application test should be conducted to ensure that the material withstands the specific requirements.

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Additional:

Instructions for use, chemical resistance lists, and quality certificates for each product are available for download from www.tpp.ch.

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Disclaimer:

TPP products are intended solely for use in research and are not suitable for clinical, diagnostic, or therapeutic applications.

TPP assumes no liability for damages or defects resulting from improper or unauthorized use. Users are responsible for storing, handling, and using the products in accordance with the provided instructions.

TPP's recommendations are general guidelines and may not cover all possible scenarios. TPP is not liable for indirect, incidental, consequential, or special damages arising from the use or misuse of these products.

Swiss law governs these terms of use and any resulting legal matters.

Literature:

- [1] Amanda Capes-Davis, R. Ian Freshney (2010) Freshney's Culture of Animal Cells: A Manual of Basic Technique and Specialized Applications (8th Ed.) - Wiley (p.180)