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Instruction for Use TPP Serological Pipette

TPP's sterile serological pipettes are for dispensing aqueous fluids such as cell culture media.

The color-coding on the pipette and dispenser allows for easy identification.

The negative graduation and reservoir provide more working volume than conventional serological pipettes.

The specially designed dispenser with perforated peel-off tab serves as a dispenser for easy collection directly at the bench.

The TPP serological pipette is for single use only.

Safety instructions:

1. Handling and Safety

When dealing with biological material, adherence to national regulations is mandatory. Furthermore, it is essential to ensure compliance with the established biological safety level of the laboratory, along with the safety data sheets and the manufacturer's instructions for use.

Wear your personal protective equipment.

2. Risk of contamination

Always follow the rules for aseptic techniques and good laboratory practices. Open the packaging only immediately before use and use only items that are visually flawless. Dispose of any products that show visible defects or signs of contamination.

3. Storage

Store the TPP products at room temperature and protect them from direct UV radiation.

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Instruction

- *Check the expiration date (EXP) on the label and packaging. Only use products with a valid EXP date.*
- *Before use, verify that the packaging is intact, as the consumable is only considered sterile if the packaging is undamaged.*
- *Open the packaging in a sterile environment and remove the product.*
- *Tear open the wrapper to expose the mouthpiece of the pipette.*
- *Place the mouthpiece into the pipetting-aid, keeping the pipette within its wrapping to avoid contamination.*
- *If the pipette does not fit properly, replace the silicone pipette holder inside the nosepiece of the pipet aid.*
- *With the pipette in place, hold the aid in place and then remove the wrapper.*
- *Use the pipetting aid to fill the pipette to the desired volume by aspirating up to the lower meniscus.*
- *Do not overfill the pipette! Do not allow liquid to come into contact with the cotton plug or getting into the pipetting aid.*
- *To ensure accuracy, hold the pipette vertically in front of the eye and measure the volume at the base of the meniscus.*
- *Allow for parallax error.*
- *Position the pipette in the target vessel and allow the liquid to flow in slowly (on the vessel wall in the case of centrifuge tubes).*
- *Before regular use of the product, perform a test run with the chosen settings to verify its suitability.*

Advice and Recommendations:

- *If pipetting in a sterile environment and the pipette touches a non-sterile surface, discard the pipette and replace it with a sterile one.*

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Additional:

Instructions for use, chemical resistance lists, and quality certificates for each product are available for download from www.tpp.ch.

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Packaging unit for triple bagged TPP Serological Pipettes

Nr. No.	Karton Case	Produktbeschreibung 3-fach verpackt Product Description Triple Bagged	Stk. / Beutel Qty. / Bag	Stk. / Karton Qty. / Case
394525	1/4	Serologische Pipette 25 ml mit Reservoir <i>Serological pipette 25 mL with reservoir</i>	1 bag of 3	36
394550	1/4	Serologische Pipette 50 ml mit Reservoir <i>Serological pipette 50 mL with reservoir</i>	1 bag of 3	36

Chemical Resistance of the Ink

+ = Good chemical resistance	± = Good to limited chemical resistance	- = Poor chemical resistance
Substance does not cause damage	Long-time and/or continuous exposure to the substance causes minor damages	Not suitable for exposure to the substance, damage may occur (deformation, discoloration, cracking, dissolution).

Solution	+ / ± / -	Solution	+ / ± / -
Acetic acid	±	Formaldehyde	±
Acetone	-	Glutaraldehyde	-
Acetonitrile	-	Glycerol	+
Ammonia solution	±	Isoamyl alcohol	±
Caustic soda	-	Isopropanol	±
Chloroform	-	Methanol	+
Dimethyl sulfoxide	+	Nitric acid	-
DMEM	+	Sulfuric acid	-
Ethanol	±	Triton X-100	+
Fetal bovine serum	+	Water	+

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Disclaimer:

TPP products are intended solely for use in research and are not suitable for clinical, diagnostic, or therapeutic applications.

TPP assumes no liability for damages or defects resulting from improper or unauthorized use. Users are responsible for storing, handling, and using the products in accordance with the provided instructions.

TPP's recommendations are general guidelines and may not cover all possible scenarios. TPP is not liable for indirect, incidental, consequential, or special damages arising from the use or misuse of these products.

Swiss law governs these terms of use and any resulting legal matters.