



Recommendation of appropriate operating conditions for conical TubeSpin® Bioreactors 15 / 50 mL

	Mammalian Cell culture		Insect Cells	Plant Cell culture
TubeSpin® bioreactor 15	CHO	HEK-293	SF-9	<i>N. tabacum</i> BY-2
Shaking speed rpm	200 – 300	200 – 300	200 – 300	200 – 300
Throw mm	50	50	50	50
Working volume mL	1.5 – 6	1.5 – 6	1.5 – 6	1.5 – 6


	Mammalian Cell culture		Insect Cells	Plant Cell culture
TubeSpin® bioreactor 50	CHO	HEK-293	SF-9	<i>N. tabacum</i> BY-2
Shaking speed rpm	180 – 250	180 – 250	180 – 250	180 – 250
Throw mm	50	50	50	50
Working volume mL	15 – 30	15 – 30	5 – 10	5 – 10

Max. rotational speed (50 mm shaking diameter) of different volumes of water in 15 mL TPP tubes

Timothy Bürgin^[1], Dr. Tibor Anderlej^[1]

Parameters:

- Different filling volumes of water (0.5 mL – 11 mL)
- Different setting angles of the tube holder (0°, 30°, 45° and 60°)
- 50 mm shaking diameter
- Measured by filming the shaking motion with high speed camera (in 25 rpm steps)






Mixing time of different volumes of water and different shaking frequencies in 15 mL and 50 mL TPP tubes

Timothy Bürgin ^[1], Dr. Tibor Anderlej^[1]

Method: Discoloration (with sodium thiosulfate) of amyllum solution colored with Lugol's iodine.



^[1] Kuhner AG, Birsfelden (Switzerland)



Methods from the Literature

Reference	Working Volume mL	Throw/Orbit mm	Shaking Speed rpm	TubeSpin® Bioreactor	Cells
Strnad (2010), Biotechnol. Prog.: Vol. 26 (3);pp 653-663	14 – 34	50	180 – 300	50	CHO
Xie (2011), Cytotechnology: Vol. 63 (4); pp 345-350	10 – 35	n.s	160 – 200	50	CHO
Stettler (2007), Biotechnol. Prog.: Vol. 23 (6); pp 1340-1346	5 – 10	50	140	50	CHO
De Jesus (2004), Biochem. Eng. Journal: Vol 17, No. 3; pp 217-223	5	25	200	50	CHO
Shen (2011), BMC Proc.: 5 (Suppl 8); P37	10 – 35	n.s	160 – 200	50	Sf-9
Xie (2011) Biotechnol. Letters Vol 33, pages897–902	10	n.s	75 – 90	50	Sf-9
Hacker (2018) Methods Mol Biol; 1850:123-131	10	50	180	50	CHO/HEK293
Blessing (2019) Mol Ther Methods Clin Dev. 13:14-26.	5 – 10	50	180	50	HEK293

Disclaimer

TPP products are for research use only and not for clinical, diagnostic or therapeutic use. All products are intended for use by trained personnel that are familiar with safe laboratory practices.

TPP assumes no responsibility for damage or defects resulting from improper or unauthorized use. It is the responsibility of the user to store, handle, and use the products in accordance with the instructions provided.

TPP does not warrant the completeness or accuracy of this TechDoc. TPP's recommendations are intended as general guidelines and may not cover all possible scenarios. TPP shall not be liable for any indirect, incidental, consequential, or special damages arising out of the use or inability to use the information in this TechDoc.

Swiss law governs these terms of use and any resulting legal matters.