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Instruction for Use TPP TubeSpin® Bioreactor

The TubeSpin® bioreactor is suitable for high-throughput screening, cultivation and optimization processes of prokaryotic (aerobic bacteria) and eukaryotic (yeasts, algae, HEK, CHO, Sf-9) suspension cells.

The TubeSpin® bioreactor is available with a filter screw cap with 5 or 10 openings over a PTFE filter membrane.

The TubeSpin® Bioreactor is intended for single use only.

Note:

Follow national regulations when handling biological material, use appropriate protective clothing

Take note of the regulations of aseptic working during the work process.

Instructions:

- *Check the packaging material and the product for intactness. Use only faultless material.*
- *Check the expiry date (EXP) on label and packaging. Only use products with valid EXP.*
- *Open the packaging in a sterile environment and remove an entire system for use.*
- *Open the bioreactor and fill it with the medium and inoculum according to your laboratory routine. Please refer to the optimal filling quantity (see technical data).*
- *Close the filter screw cap with ventilation holes. A continuous gas exchange takes place through the integrated hydrophobic membrane 0.22 µm.*

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- *If the aeration/gas exchange needs to be reduced, individual openings may be closed using lab appropriate labels or tape (picture 1)*

Picture 1

- *Cultivate: Shake cells in appropriate incubation shakers. Refer to the cultivation conditions of the cells (temperature, CO₂ and humidity control, shaking frequency)*
- **Shakers:** *Follow the safety instructions and operating instructions of the shaker manufacturer.*
- **Centrifuging:** *The use of suitable rotors or centrifuge adapters is recommended. Follow the safety instructions of the centrifuge manufacturer.*
- *The mechanical strength of the tubes during centrifugation is influenced by*
 - *Shape and material*
 - *Accuracy of fit in the centrifuge adapter*
 - *Temperature, centrifugation time, g-number*
 - *Chemical and physical properties of the centrifuged material*
 - *Rotor type: fixed-angle or swing-out rotor*
- *Test the multitude of influencing factors under routine conditions beforehand.*
- *Storage before use at room temperature, protected from UV light*

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Remarks and Recommendations:

- *The graduation serves as a guide for filling quantity and is no absolute value.*
- *Alcohol 90 % with mechanical stress (rubbing) might dissolve or blur the white color of the print.*
- *For freezing TPP recommends transferring samples to the TPP cryo tubes. If TubeSpin® bioreactor storage in minus temperature range is desired, TPP recommends testing the product beforehand according to your conditions.*
- *TPP assumes no liability for damage caused by, during or due to the use of centrifuge tubes in freezing processes.*

Technical Data:

Materials

Screw Cap	PE
Membrane	PTFE
Tube	PP

Additional:

Instructions for use, chemical resistance lists, and quality certificates of the respective products are available for download on the homepage www.tpp.ch .

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