



## Hettich GmbH & Co KG Centrifuges and Rotor Compatibility TubeSpin® Bioreactor 450

The following rotors of the centrifuge from Hettich are compatible for the TubeSpin® Bioreactor 450.



**ROTANTA 460**



**ROTANTA 460 R**

<b>Swing-Out Rotor</b>			<b>5699-R</b>
Hanger with cup cup bio-safe 4883	Hanger without closure	Adapter	Capacity
4890	4895	4449	4 x 87450

For further technical information on centrifuges, rotors and adapters from Andreas Hettich GmbH & Co KG please visit <https://www.hettichlab.com> and/or contact [info@hettichlab.com](mailto:info@hettichlab.com)



**Technical Data**  
**TubeSpin® Bioreaktor 450 / 600**

<b>Measurements</b>	<b>87450</b>	<b>87600</b>
<i>Volume graduation mL</i>	450	600
<i>Length mm</i>	145	183
<i>Diameter mm</i>	96	98
<i>Max. RCF x g</i>	3'500	3'500
<i>Form</i>	<i>konisch / conical</i>	<i>konisch / conical</i>
<i>Optimal filling volume mL</i>	>250	>400
Shaker:		
Recom. orbit/shaking throw mm	50	50
Shaker: Recom. speed rpm	150	150
Screw cap „solid“	87355	87355
Screw cap „filter“	87356	87356

**Disclaimer**

TPP products are for research use only and not for clinical, diagnostic or therapeutic use. All products are intended for use by trained personnel that are familiar with safe laboratory practices.

TPP assumes no responsibility for damage or defects resulting from improper or unauthorized use. It is the responsibility of the user to store, handle, and use the products in accordance with the instructions provided.

TPP does not warrant the completeness or accuracy of this TechDoc. TPP's recommendations are intended as general guidelines and may not cover all possible scenarios. TPP shall not be liable for any indirect, incidental, consequential, or special damages arising out of the use or inability to use the information in this TechDoc.

Swiss law governs these terms of use and any resulting legal matters.