




Mounting & Disassembling a TPP Tube Rack



Rack "S" and "L" for centrifuge tubes

- Identification: Alphanumeric coded
- Storage: Space saving, as it can be quickly disassembled
- Holds:
 - Rack "S": 18 x 15 mL and 10 x 50 mL
 - Rack "L": 30 x 15 mL and 20 x 50 mL

Measurements	99017	99019
Length mm	209	209
Width mm	98	172
Height mm	62	62
Tubes per rack	18	30
Materials	PP	
Working Temperature	+121 to -190 °C	
Autoclavable	yes	
Sterilizing instructions:	Please sterilize disassembled only	



Disclaimer

TPP products are for research use only and not for clinical, diagnostic or therapeutic use. All products are intended for use by trained personnel that are familiar with safe laboratory practices.

TPP assumes no responsibility for damage or defects resulting from improper or unauthorized use. It is the responsibility of the user to store, handle, and use the products in accordance with the instructions provided.

TPP does not warrant the completeness or accuracy of this TechDoc. TPP's recommendations are intended as general guidelines and may not cover all possible scenarios. TPP shall not be liable for any indirect, incidental, consequential, or special damages arising out of the use or inability to use the information in this TechDoc.

Swiss law governs these terms of use and any resulting legal matters.